# RapID™ SS/u System

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## RapID™ SS/u

For the rapid ID of commonly isolated urinary tract microorganisms from human specimens. Panel includes 11 substrates for the identification of 12 organisms.

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### RapID™ Systems and Accessory Products

All systems are packaged 20 panels/pk. Color guides for each system are available upon request from your Technical Sales Representative or our Technical Service Department.

Accessory Products	Packaging	REF
RapID™ Inoculation Fluid	1ml x 20 tubes/Pk	R8325102
RapID™ Inoculation Fluid	2ml x 20 tubes/Pk	R8325108
RapID™ Spot Indole Reagent	15ml/Btl	R8309002
RapID™ Nitrate A Reagent	15ml/Btl	R8309003
RapID™ Nitrate B Reagent	15ml/Btl	R8309004
BactiDrop™ Oxidase	0.75ml x 50 ampules/Pk	R21540
McFarland Equivalence Turbidity Standard No. 1	Each	R20411
McFarland Equivalence Turbidity Standard No. 2	Each	R20412
McFarland Equivalence Turbidity Standard No. 3	Each	R20413
McFarland Equivalence Turbidity Standard No. 4	Each	R20414
ERIC® (Electronic RapID Compendium)	Each	R8323600

RapID™ System	Inoculation Fluid <sup>1</sup>	McFarland Standard No. <sup>2</sup>	Spot Indole	Nitrate A	Nitrate B	Oxidase	Culti-Loop® QC Set <sup>3</sup>	Paper Code Compendium	Available on ERIC®
ANA R8311002	1ml	3	√	–	–	–	R4653050	R8322002	√
CB PLUS R8311008	2ml	4	–	√	√	–	R4653048	R8328008	√
NH R8311001	1ml	3	√	√	√	√	R4653051	R8321001	√
NF PLUS R8311005	1ml	1	√	√	–	√	R4653054	R8325005	√
ONE R8311006	2ml	2	√	–	–	√	R4653056	R8328006	√
STAPH PLUS R8311009	2ml	3	–	√	√	–	R4653047	N/A	√
SS/u R8311004	1ml	1	√	–	–	–	R4653053	R8324004	√
STR R8311003	1ml	1	–	–	–	–	R4653052	R8323003	√
YEAST PLUS R8311007	2ml	N/A	–	–	–	–	R4653060	R8327007	√

1. RapID Inoculation Fluid is formulated for preparation of bacterial inoculums, which are adjusted to specific densities using McFarland Equivalence Turbidity Standards and then inoculated into an appropriate RapID Identification System panel.

2. McFarland Equivalence Turbidity Standards are accurate, stable latex turbidity standards recommended for use with RapID systems. These convenient turbidity standards are packaged in a plastic case with a visual comparison card.

3. Culti-Loop® QC sets and kits are recommended for consistent, convenient quality control of the RapID systems. Culti-Loops are ready-to-use, stabilized, preserved, viable microorganisms in a disposable inoculating loop. The organisms can be rehydrated in liquid or streaked directly onto an appropriate eger medium.

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### ERIC®

ERIC® (electronic code compendium) is a complementary computerized database to support all RapID Systems.

- Windows®-based software system
  - Computerized code compendium to process all RapID system microcodes
- User friendly format for RapID code interpretation
- Expanded databases with additional organisms and updated nomenclature
  - Increased identification capabilities
  - Reliable, comprehensive results

ERIC includes all the compendium information needed including: probability percent, bioscores, contraindicated test results, nomenclature, and clinical relevance comments. It also provides save and report option to track identification patterns and organism trends.

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### RapID™ Color Guide

Results of each RapID test can be compared to a convenient RapID Color Guide. The easy-to-read guides provide reference information at a glance. Complete instructions for use are included in every RapID System.

#### RapID™ ANA II Color Guide

992-143

Test	Cavity	Positive Reactions	Negative Reactions		
URE	1				
BLTS	2				
aARA	3				
ONPG	4				
aGLU	5				
BGLU	6				
aGAL	7				
aFUC	8				
NAG	9				
PO4	10				
LDY	3				
GLY	4				
PRD	5				
PAL	6				
ARG	7				
SER	8				
PYR	9				
IND	10		Any Other Color		

#### Comments

Without the addition of any reagents, read cavities 1-10 and record results.

**Cavity 1:** Development of a red or purple color is a positive test, a yellow or orange color is a negative test.

**Cavities 2-10:** Development of a medium to bright yellow color is a positive test, a clear, tan, or pale yellow color is a negative test.

**Add RapID ANA II Reagent to cavities 3-6. Allow at least 30 seconds but no longer than 2 minutes for color development.**

**Cavities 3-9:** Development of a distinctive purple, red, or dark pink color is a positive test, a pale color is a negative test.

**Add Spot Indole Reagent to cavity 10.**

**Cavity 10:** Development of a blue or blue-green color is a positive test, any other color is a negative test.

Notes: The RapID™ Color Guides are provided as an educational aid to be used in conjunction with the instructions for use for the product. The results shown above are the only approved color swatch combinations for positive and negative results.

Min. 2014 © Remel  
992-143 (2014) (2015) (2016)  
992-143-01-00 (2014) (2015) (2016)

14386 Series 10-10-14  
Lot# 14386-10-10-14

**remel**  
Remel Science Products

## RapID™ Systems plus ERIC® equals results

Remel RapID™ systems are ideal for manual testing and companion confirmation tests for automated systems.

- Accurate identification of more than 400 medically important organisms in 4 hours.
- Common, simple, one-step inoculation with aerobic incubation.
  - Less than 3 minutes required to set up and read results.
  - No oil-overlays necessary.
  - Compatible with a wide range of media formulations and culture types.
- Test results based on easy-to-read color reactions in each individual well.
  - Color-based enzymatic cleavage of specific chromogenic and conventional substrates for accurate results.
  - Convenient RapID provide interpretation at a glance.
  - ERIC® software is reliable, comprehensive, and easy to use..
- Complete set of accessory products designed specifically for use with RapID Systems for optimum performance.

## Product Comparison Guide

	RapID™	Others
Incubation time	4 hours	Range from 4 - 48 hours
Inoculation procedure	1 step	Multiple pipette steps
Aerobic incubation	All	Some
Windows-based compendium	Yes	Some