

Certificate of Analysis

Product Code: ETO4/MD
Product Name: OXIDASE STRIPS GERMAN
Lot Number: 503201
Expiry: 30/06/2024

Quality Assurance

Mast products are manufactured in accordance with an ISO9001/13485 Quality system which covers design, manufacture, and supply of in-vitro diagnostic devices and consumables.

Quality Control Testing

The above products were tested in the Quality Control Laboratory and met defined acceptance criteria using our standard test method.

Performance: Tests for bioburden, uniformity and microbiological performance were within specification.

SATISFACTORY FOR ISSUE



Susanne Chapman
Quality Control Laboratory Manager