

Certificate of Analysis

Product: CombiScreen® Drop Check **Art.-No.:** 93015
Filling: 2 x 5 ml

Lot number: A 130

Manufacturing date: 2024-07-01

Expiry date: 2026-01-31

Storage: Store between +2 and +8 °C.

Item	Specifications	Results	QC
Target Ranges see package insert.			
Solution test			
Level 1	Slightly yellow solution	slightly yellow	passed
	No particles or sediment (clear solution)	clear	passed
Level 2	Slightly brown solution	slightly brown	passed
	No particles or sediment (clear solution)	clear	passed
Material test			
Tube cap	Cap closes tight	tight	passed

This lot complies with our specification.

Date of release: 2024-08-20

This certificate of analysis has been automatically generated and is valid without signature.

Elena Schwietz
Manager Quality Control