

## Certificate of Analysis

**Product:** CombiScreen® Dip Check **Art.-No.:** 93010  
**Filling:** 2 x 15 ml  
**Lot number:** A138  
**Manufacturing date:** 2025-09-12  
**Expiry date:** 2027-04-30  
**Storage:** Store between +2 and +8 °C.

| Item                              | Specifications                               | Results         | QC     |
|-----------------------------------|----------------------------------------------|-----------------|--------|
| Target Ranges see package insert. |                                              |                 |        |
| <b>Solution test</b>              |                                              |                 |        |
| Level 1                           | Slightly yellow solution                     | slightly yellow | passed |
|                                   | No particles or sediment<br>(clear solution) | clear           | passed |
| Level 2                           | Slightly brown solution                      | slightly brown  | passed |
|                                   | No particles or sediment<br>(clear solution) | clear           | passed |
| <b>Material test</b>              |                                              |                 |        |
| Tube cap                          | Cap closes tight                             | tight           | passed |

This lot complies with our specification.

**Date of release:** 2025-10-08

This certificate of analysis has been automatically generated and is valid without signature.

Elena Schwietz  
Manager Quality Control